

Glass vials for pharmaceutical applications

We offer a wide variety of vials from 1 to 50 ml in different designs, with or without blowback (European and American version) which covers any need our customers may have. Production according to industry standards or your specifications. Packaging is performed under controlled environmental conditions.

Our product range includes clear and amber vials in glass types I, II, III.
Our vial assortment comprises:

- **Pharmaceutical vials**

- Injection vials (serum vials)
- Lyophilisation vials
- Tablet vials
- Large OD vials

- **Diagnostic vials**

- Screw thread vials
- Screw thread tubes
- Chromatography vials

- **Special vials**

- ADD-Vantage® vials
- Two-compartment vials
- Sampler vials
- Onion skin vials

We offer the following options enhancing glass quality and processability:

- Ammonium sulphate treatment (AST, see page 6.2)
- Optimised designs for lyophilisation
- High quality pharmaceutical serum vials using proprietary visual inspection systems
- Low alkalinity option
- Siliconised vials
- Our Pharma Plus line to meet the most stringent cosmetic and dimensional control requirements

